
EDOS® 5221

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Eppendorf -
Netheler - Hinz GmbH
22331 Hamburg
Tel.: (0 40) 5 38 01-0
Teletex: 40 30 61 = EGHAM
Fax: (0 40) 5 38 01 556

5221 900.105-00/1092/Mx
Printed in Germany

eppendorf

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The screenshot shows a device display with the following elements:

- Top bar: Disp Pip Mix Titr Dil Manu ◊Pgm ◊Opt
- Tip: 2.50 m
- Navigation: >> MULTI DISPENSE <<
- Index bar: 1 2 3 4 5 6 (with '1' highlighted)
- Speed controls: Speed ↑ (3 bars) and Speed ↓ (3 bars)
- Disp. (1): 50 µl
- Rep. (1): 3 / 2
- Sum: 2050 µl

Annotations on the right side of the image:

- A horizontal line with an asterisk (*) points to the Index bar.
- A vertical line with an exclamation mark (!) points to the Rep. (1) value.

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Liquid:	Bidistilled water.
Reference temperature:	20 – 25 °C, constant to ± 0.5 °C.
Number of determinations:	15, in accordance with DIN 12 650 with Original Eppendorf pipette tips.
Main power connection:	230 V/50 Hz
Power requirement:	20 W
Fuse for main power supply:	primary for 230 V 0.2 A time-lag secondary: 1.0 A time-lag

Technical specifications subject to change!

Delivery package:

- 1 control unit
- 1 dispensing grip
- 1 power supply
- 2 pipetting adapters (Pip 100 + Pip 1000)
- 1 rest
- 1 Combitip in each size
- 1 adapter for 25 and 50 ml Combitips
- 1 main power cable
- 1 set of fuses
- 1 instruction manual

11 Technical data

EDOS meets the requirements set by DIN 12 650 with all relevant dispensing components.

For Pip 10, 100, 1000, Multi 50 and Multi 200, the data applies to pipetted volumes; for Combitips, it applies to dispensed volumes.

When dispensing in free jet with Combitip ≤ 2.5 ml, a 100 μ l (10 μ l) pipette tip must be attached. With Combitips ≤ 1.25 ml, the pipette tip must be held against the inside of the vessel.

Dispensing component	Inaccuracy %			Imprecision %		
	1/50 vol.	1/10 vol.	Max. vol.	1/50 vol.	1/10 vol.	Max. vol.
Pip 10		± 2.5	± 1.0	1.8	0.4	
Pip 100		± 1.0	± 0.6	0.4	0.2	
Pip 1000		± 0.8	± 0.6	0.3	0.2	
Multi 50	in preparation			in preparation		
Multi 200	in preparation			in preparation		
0.05 ml Combitip	± 2.5	± 2.5	± 1.0	8.0	2.5	1.0
0.5 ml Combitip	± 0.9	± 0.9	± 0.9	1.2	0.7	0.3
1.25 ml Combitip	± 0.9	± 0.9	± 0.5	0.8	0.5	0.2
2.5 ml Combitip	± 0.8	± 0.8	± 0.5	0.8	0.5	0.2
5.0 ml Combitip	± 0.8	± 0.8	± 0.5	0.6	0.3	0.15
12.5 ml Combitip	± 0.5	± 0.5	± 0.4	0.4	0.25	0.15
25.0 ml Combitip	± 0.4	± 0.3	± 0.3	0.6	0.35	0.15
50.0 ml Combitip	± 0.3	± 0.3	± 0.3	0.35	0.25	0.15

The **piston sealing rings** are made of chemically resistant fluoropolymer, the piston of chemically inert ceramics.

The seal can be exchanged separately when **Pip 1000** is used.

The seal of the **Pip 100** is located in the lower housing of the pipette, the seal of the **Pip 10** in the nose cone. Both seals can only be exchanged together with the adapters (Order no. for lower housing of the pipette for Pip 100: 5221 501.006).

In addition, all parts marked in the figure by an order number are available as spare parts.

When the piston has been cleaned with alcohol, it must be lubricated slightly (silicone grease, Sec. 12, see Ordering information).

If you attempt to pull a dispensing component out of the dispensing grip without moving the piston mounting to the Zero Position, the piston can slide out of the mounting. Reinsert the piston as follows:

Set the lowest dispensing speed (1 block).

Carry out a Reset. If the motor makes a noise at the end of the Reset (groans), stop by re-pressing the R key on the dispensing grip. Push the piston back into the mounting by rotating slightly.

Important note:

Please only use the accessories recommended by Eppendorf. Using spare parts or disposables which we have not recommended can reduce the precision, accuracy and life of our instruments. We do not honor any warranty or accept any responsibility for damage resulting from such action.

10 Maintenance, cleaning, spare parts


Fig. 11:

- 1 = Ejection sleeve
- 2 = Nose cone
- 3 = Piston seal
- 4 = Sealing spring
- 5 = Piston
- 6 = Lifting rod
- 7 = Clutch
- 8 = Lower housing of pipette

The device is maintenance-free. Use diluted rinsing or soap solutions to clean all plastic parts. Disinfectants with an alcohol base must not effect the plexiglass cover of the display in the relative long term. Wipe cover carefully after decontamination.

All pipetting adapters are completely autoclavable.

The pipetting adapters can be taken apart (see Fig. 11):

Pip 1000: Unscrew clutch. Pull out piston and lifting rod. Remove piston seal from the piston. Pull lifting rod off the piston.
Unscrew nose cone (normal screw thread). Pull off ejection sleeve.

Pip 100 + Pip 10: Unscrew clutch. Pull out piston and lifting rod. Pull lifting rod off the piston.
Pull off ejection sleeve. Unscrew nose cone (normal thread).

Error message/ Error	Cause	Remedy
Dispensing grip makes noise "groans" when ▽ pressed.	- Motor moves towards Zero Position due to incorrect step counting (dispensing component sluggish).	Stop motor by pressing R. Remove dispensing component. Execute Reset. Clean or exchange dispensing component.
Error message "Eject manually".	- Pipette tip fitted too tightly on pipetting adapter.	Eject tip manually. (Caution: Do not pull pipetting adapter out to side!)
Error messages "VOLUME-ERROR" "STEPPER-ERROR"	- Combipip jammed because of attempt to remove it when EDOS not in home position.	Press in Combipip. Press Reset.
	- Speed too high for viscous solution.	Reduce speed.
	- Electronics error.	*Call service.
Error messages "GRIP KEY-ERROR" "GRIP EE-ERROR" "GRIP CAL-ERROR" "EEPROM-ERROR" Continuous beeping.	- Electronics error.	*Call service.

*** General note:**

If error occurs, switch off EDOS to reset error message. If error re-occurs after switch-on, call service. Disconnect EDOS power unit from the power supply.

9 Troubleshooting

Error message/ Error	Cause	Remedy
Dispensing component cannot be inserted.	- Piston of the dispensing component not pushed in completely.	Push in piston completely.
	- Piston mounting not in Zero Position.	Execute Reset. Hold down R key when inserting pipetting adapters.
Combitip drips.	- Combitip leaks.	Exchange Combitip.
	- Attached tip is loose.	Press on tightly.
Pipette tip drips.	- Tip fits loosely on the nose cone of the pipetting adapter.	Press on tightly.
	- Wrong pipette tip.	Use Original Eppendorf pipette tip.
	- Nose cone loose.	Tighten nose cone.
	- Nose cone, O-ring of nose cones (Multi 50 and Multi 200), piston sealing-ring, piston damaged.	Exchange defective parts.
Error message "Fill!".	- ▽ pressed although not sufficient liquid for dispensing.	Press Δ to fill again.
Error message "Is full!".	- Δ pressed, although Combitip/pipette tip is full.	Press ▽ to empty.
Error message "Press down!"	- Δ pressed, although blow-out must be performed.	Perform blow-out.
Error message "Use another tip!"	- Function cannot be performed with the dispensing component selected.	Insert suitable dispensing component.
Error message "Remove tip!"	- The control unit can distinguish between pipetting adapters and Combitips using a sensor built in the dispensing grip. This error message is displayed when the wrong dispensing component is used.	Insert correct dispensing component (Combitip or pipetting adapter).

8 Options

Options (*Opt* function)

After activation of the *Opt* function, you can choose between:

Beep: For switching the beeper off and on (with Up/Down cursor keys).

Key click: If Key click is selected, every keypress at the control unit and the reaching of the target position is confirmed by an acoustical tone.

Calibr.: With this calibration factor, the piston stroke for liquids whose viscosity and density differ from that of water is corrected. For this purpose, determine the volume of the liquid dispensed gravimetrically: divide weight of the liquid by its density. Divide the expected nominal value of the liquid volume by the volume value calculated and use this ratio as a calibration factor.

The display "Cal : #####" appears in the entry field when a calibration factor other than 1.000 has been entered.

Remote: The serial interface for remote control of the device can be set here.

Exit: By selecting and confirming Exit, access to the next higher display level in the command line is attained.

7.2 Selecting a set of parameter values

Note: Until the final selection of a set of parameter values (by point 7) **Enter**, see below, the set of parameter values last entered or selected is retained (activated). Even if, after points 4), 5) and 6), other sets of parameter values are shown in the display, the last set of parameter values entered/selected appears in the display after interruption of the **Recall** function by pressing Mode.

By this means, it is possible to select a suitable, overwriteable memory location before storing a set.

Key sequence at control unit	Meaning
1) Mode	Activation of the command line.
2) ◀▶ cursor keys	Marking of Pgm .
3) Enter	Activation of the subcommand line Pgm. Recall is marked.
4) Enter	Activation of the memory.

The field for the entry of the memory location number is marked in the bottom line of the display (Source Nr.: #).

5) ◀▶ cursor keys	Choice of function.
6) ▲▼ cursor keys	Choice of memory location number.
7) Enter	Call-up of the set of parameter values displayed in the entry field. 2-second comment "Nr. # recalled!"

After point 7), the control unit moves back to the function mode and the **function selected** can be **executed** immediately with the dispensing grip.

The set of parameter values activated for the last function is displayed in the entry field **until point 4)** is carried out.

After point 4), the set of parameter values last stored or selected **for this function** with **Pgm** is displayed in the entry field.

After point 5), the set of parameter values of a new function selected with the right/left cursor keys is displayed in the entry field, i.e. the last set stored/selected with **Pgm**.

After point 6), the set of parameter values belonging to the memory location number of this function is displayed in the entry field.

or interruption:

instead of points 4), 5) or 6), press "Mode":

Comment: "Cancelled" (2 seconds long), return to point 3). After interruption, the set of parameter values last activated appears immediately in the entry field.

From point 3) return to point 1) by re-pressing "Mode".

7 Storing and selecting the programs

Storing and selecting the programs (*Pgm* function)

All functions *last* defined for the different functions and dispensing components at the entry field are stored after switching off the control unit.

It is also possible to file up to 30 parameter sets *per function* in a program memory. A calibration factor is stored for every set as an additional parameter (see Sec. 8).

7.1 Storing a set of parameter values

After a set of parameter values for a chosen function have been entered (see Sec. 3.2.1 and 3.2.2), they can be stored:

Key sequence at control unit	Meaning
1) Mode	Activation of the command line.
2) ◀▶ cursor keys	Marking of Pgm .
3) Enter	Activation of the subcommand line Pgm .
4) ◀▶ cursor keys	Marking of Store .
5) Enter	Activation of the memory.

The field for the entry of the memory location number is marked in the bottom line of the display (Dest. Nr.: #).

6) ▲▼ cursor keys	Entry of the memory location number
7) Enter	Storage of the set of parameter values under the memory location number. 2-second comment "Nr. # stored!" Return to point 4).

or

7) Mode	Interrupted, return to point 4). 2-second comment "Cancelled!"
8) Mode	Return to point 1).

During execution of the **Store** function, the set of parameter values last entered/selected which should be stored is displayed in the entry field.

6.6.2 Acoustical signals:

One long tone marks the end of a dispensing sequence.

Two consecutive tones serve as a warning signal (e.g. wrong key pressed).

The acoustical signals can be switched off (see Sec. 8).

Every keypress on the control unit and the dispensing grip is confirmed by a tone when "Key click" is switched on (see Sec. 8).

6.6.3 Comments

The status of the device is described by comments which appear in the bottom line of the display field on the left:

Finished!:	A sequence with preselected repeat number is finished.
Multi finished!:	A Multi Dispense sequence is finished.
Zero Position:	The piston is pushed completely into the dispensing component.
Resetting...:	Reset is executed.
Blowing out...:	Blow-out is executed.
Stopped:	After stopping the piston by pressing R or Δ .
Select new item!:	The command line is activated.
Press Reset!:	After switching on the device, entry of a new dispensing component, stopping piston movement by pressing R.

Additional notes: see Sec. 9.

6.6 User prompting

6.6.1 Optical signals

The keys of the dispensing grip light up to indicate the next key to be pressed.

The Δ key is **green** and flashes when refilling must take place for the processing of a dispensing sequence already begun.

A key *must* be **pressed** when:

- the R key is **red**,
- the ∇ key is **yellow** because a blow-out is necessary (**Pip**, **Mix**, **Dil** functions).

A key *can* be **pressed** but does not have to be when:

- Δ key **green** and flashes Piston position of a pipetting adapter = **Pip home position** (after execution of blow-out). Instead of filling, blow-out can be performed as often as desired by pressing ∇ .
- ∇ key **yellow** Indicates that liquid can be dispensed. The Combitip or tip is not completely full (**Disp**, **Tit** and **Manu** functions).
A Combitip can be filled with the liquid volume for several dispensing sequences (by pressing the Δ key several times).
The end of every sequence from a Combitip thus filled is indicated by the **green** light of the Δ key and a long tone.

6.4 *Speed of piston movement* ("Speed parameter")

The setting of a medium speed (4 blocks) means:

- for pipetting adapters, a speed of liquid movement which is the same as during normal handling of manual pipettes,
 - for Combitips, liquid movement which enables dispensing of aqueous solutions with exact drip separation in free jet.
- Exceptions: see Sec. 11.

When pipetting viscous liquids, lower speeds may be appropriate. With the *Mix* function, higher ones may be appropriate.

Reset is carried out at the speed set for dispensing.

Tip ejection after pressing the R key (pipetting with pipetting adapters) is carried out at a defined, fixed speed.

6.5 *Entering parameter values during execution of a function, new start of a repeat counting*

During execution of a function, the values for the *Speed* parameter can be changed as soon as the piston stops moving.

New values for the *Volume* and *Repeats* parameters can be entered at any time during the *Disp* function when the piston is at a standstill, regardless of the position of the piston.

A new value can only be entered for the *Volume* and *Repeats* parameters with the *Pip*, *Mix* and *Dil* functions, if the piston of the dispensing component is in the *Pip* or *Disp* home position (see Sec. 6.1 and 6.2).

The *setting of another dispensing component* results in:

- the display of the set of parameter values last entered for this dispensing component,
- the prompt to Reset.

The *Repeats count* always begins:

- after a change in a volume value or a repeat number,
- after "Mode" and "Enter" are pressed in succession.

6.3 Controlling the piston movement, Reset

6.3.1 Stopping piston movement

By pressing the R key *briefly*, the piston can be stopped during the upward or downward movement (emergency stop function). With the *Disp*, *Plp*, *Mix*, *Dil* functions, it is then necessary to perform a Reset.

With the *Disp* function, filling can be discontinued by pressing the Δ key of the dispensing grip again (and continued by re-pressing).

A discontinued dispensing step is not counted. As long as a Combitip contains enough liquid for the dispensing of a partial volume, this is dispensed and counted.

Note: With the *Disp* function, *no* relaxation stroke is executed when filling is discontinued by pressing R briefly.

6.3.2 Piston movement with the *Tit* and *Manu* functions

With the *Tit* and *Manu* functions, the piston only moves as long as the $\nabla\Delta$ keys are held down.

6.3.3 Reset function

Press the R key for a while to initiate an immediate Reset. Press the R key again to bring the piston to a standstill. Once a Reset is started, it must be completed.

Reset has the following functions:

- to move down the piston mounting so that dispensing components can be inserted,
- to discontinue a dispensing process,
- to define the Zero Point of the dispensing spindle in the dispensing grip by orienting on a calibration point.
By this means, the accuracy of dispensing is assured.

Except after interruption of the *Disp*, *Plp*, *Mix* and *Dil* functions by pressing R, a Reset is required when:

- the device is switched on,
- a dispensing component is exchanged.

This is prevented by the **residual stroke lock** of the Multipette and EDOS. The last 1/52 of the Combipip volume is not utilized. The piston of the Combipip is not moved deeper than the **Disp home position** when ∇ is pressed (Fig. 10b).

The Combipip is filled from the **Disp home position** to continue dispensing (Fig. 10c).

After filling, the piston always automatically executes a so-called **relaxation stroke** of 1 mm = 1/52 of the Combipip length: the piston moves past the point calculated for filling by this amount and then back to the point. This is necessary to bring the piston mounting into contact with the dispensing spindle before the down-movement. The relaxation stroke corresponds to the first dispensing step to be discarded when dispensing with the Multipette.

If a 100 μ l pipette tip must be attached to the Combipip to achieve the technical data, proceed according to the above description. This also applies if the Combipip/pipette tip is held against the inside of the vessel during dispensing (see Sec. 11).

6.2 Liquid transfer with pipetting adapters Pip 1000, Pip 100, Pip 10, Multi 200, Multi 50.

When using these dispensing components, proceed as for manual pipettes with blow-out, see Sec. 5.2.

Part of the stroke distance must be "reserved" for **blow-out** and **tip ejection** before filling.

The piston mounting therefore moves automatically from the Zero Position to the **Pip home position** when the R key of the dispensing grip is released.

From the Pip home position, the blow-out position is reached by pressing the ∇ key.

The blow-out is initiated and executed with EDOS in the same way as with a manual pipette: when ∇ on the dispensing grip is pressed, the piston of the pipetting adapter moves down to the blow-out position and does not move back until the ∇ key is released.

A **relaxation stroke** is performed after every upward movement of the piston when pipetting adapters are used and the **Manu** function is selected (or initiated).

With the **Pip**, **Mix** and **Dil** functions, the piston moves a little below the Pip home position during dispensing so that correction is performed. The piston then moves to the Pip home position during blow-out.

A **blow-out** must be performed with the **Pip**, **Mix** and **Dil** functions.

6 Useful hints


Fig. 10: Liquid transfer with a Combitip

a = Aspirating beginning with the Zero Position (Z)

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* = Air bubble

Read the following information to intensify your knowledge on the functions and options of EDOS.

6.1 Liquid transfer with a Combitip

Using a Combitip and the Eppendorf Multipette, you can dispense according to a method which is tried and tested throughout the world (see Fig. 10):

Liquid is aspirated into the Combitip from the Zero Position (piston of the Combitip pushed completely into the cylinder after Reset, Fig. 10a).

After aspirating, a small air bubble is **always** under the piston base because the contour of the piston base never completely fills out the contour of the cylinder base.

This air bubble does not disturb the precision of dispensing.

However, it would cause an error if the piston is pushed down to the base of the cylinder during dispensing and an aliquot of air is dispensed instead of liquid.

5.6 Aspirating and dispensing any amounts of liquid (*Manu* function, manual mode)

This function can be executed with all dispensing components.

Volume displays at the control unit:

Piston position **SUM**: Volume **aspirated** = for Combipips, volume between piston base and Zero Position, see Sec. 4.1 + 6.1; for pipetting adapters, volume between piston base and **Pip home position**, see Sec. 6.2.

Rel.Vol. parameter: **Dispensable** volume (when using Combipips < **SUM** because of remaining stroke, see Sec. 6.1)

Programming

Activate **command line**, select **Manu** function.

Activate **entry** field (see Sec. 3.2.1).

Where appropriate, enter a new **Tip** and **Speed** (select parameter with the "Enter" key, set values with the cursor keys).

Executing with the dispensing grip:

Keys light up:

- 1) R key red: Press R, hold down if Pip 1000, Pip 100, Pip 10, Multi 200, Multi 50 is used. Insert dispensing component and attach pipette tip(s), where appropriate.
- 2) Δ key green: Filling: press Δ and hold down as long as liquid is aspirated.
- 3) ∇ key yellow: Liquid can be dispensed as long as ∇ is pressed. The lowest piston position which can be reached with ∇ key is for Combipips the **Disp** home position or for pipette adapters the **Pip** home position or the **blow-out position** (see Sec. 6.1 and 6.2). The Zero Position is only reached by pressing R.

Press "Mode" and then "Enter", to set the **Rel.Vol.** parameter to 0. During subsequent dispensing, the volume **dispensed** is displayed, marked by a minus sign.

Note:

When Combipips are used, the **Disp** function can be selected from the **Manu** function at any position of the piston without a Reset.

This enables the aspiration of an unknown volume of liquid into a Combipip, its measurement and dispensing (aspirated volume - remaining stroke volume) in as many steps as desired (remaining stroke: see Sec. 6.1).

Example 1:

Rinsing 10 µl sample into 2,000 µl solution, repeat 6 times.
6 x 2,000 µl solution are required.

Key sequence at control unit:

Enter Pip 100 for **Mix** function

Mix.Vol. = 15 µl

Repeats = 3

then:

Activate **command line**. Select **Dil** function.

Activate **entry field** (see Sec. 3.2.1).

After selection of **Dil** and entry of the necessary dispensing component (Pip 100):

- | | Meaning: |
|---------------------------|---|
| 1) Enter: | Confirm Pip 100.
Activate Dil.Vol. parameter. |
| 2) Cursor keys: | Set 10 µl. |
| 3) Enter: | Activate Repeats parameter. |
| 4) Cursor keys: | Set Repeats = 6 |
| 5) Enter: | Activate Speed parameter. |
| 6) ▲▼ cursor keys: | Set low speeds. |

Subsequent mixing of all dilutions, e. g.:

Mix function, Pip 1000, Mix.Vol. = 500 µl, Repeats = 3,

Speed up/Speed down = 4 blocks

Executing with dispensing grip: see above

Example 2:

Producing a dilution series from 1 : 10 to 1 : 10⁶.
100 µl should be diluted in 900 µl solution.
6 x 900 µl solution are required.

Key sequence at control unit:

Enter Pip 1000 for **Mix** function

Mix.Vol. = 500 µl

Repeats = 3

then:

Activate **command line**, select **Dil** function.

Activate **entry field** (see Sec. 3.2.1).

Enter Pip 1000

Dil.Vol. = 100 µl

Repeats = 6

Speed up/Speed down = 4 blocks

Executing with dispensing grip: see above.

5.5 Diluting (*Dil* function)

Transferring a volume into a diluting liquid and then mixing. This mixing is executed according to the set of parameter values defined under the *Mix* function.

First, define *Mix* function.

The *Dil* function can only be performed with Pip 1000, Pip 100, Pip 10, Multi 200 and Multi 50.

<i>Dil.Vol.</i> parameter:	Display of the volume to be transferred
<i>Repeats</i> parameter:	Number of dilution steps (one dilution step = transfer of volume and subsequent mixing)

Programming:

Activate **command line**, select *Dil* function.

Activate **entry field** see Sec. 3.2.1).

Enter new **Tip**, where appropriate, as well as the ***Dil.Vol.***, ***Repeats*** (0 – 99), ***Speed*** parameters (select parameter with the "Enter" key, set values with the cursor keys).

Executing with the dispensing grip:

Keys light up:

- | | |
|--------------------------|--|
| 1) R key red : | Press R key and hold down. Insert dispensing component and release R. Attach pipette tip(s). |
| 2) Δ key green : | Aspirating: press Δ key. |
| 3) ∇ key yellow : | Dispensing and mixing:
Hold tip in the liquid to be pipetted and press ∇. The sample is pipetted , all of the liquid is mixed thoroughly if suitable values are selected for the <i>Mix.Vol.</i> parameter in the <i>Mix</i> function. |

Execute blow-out (dispensing of remaining liquid from the tip):

- | | |
|---------------------------------|---|
| 4) ∇ key yellow : | Blow-out: hold tip against the inside of the vessel.
Press ∇ and hold down. Piston moves down.
Pull tip away from remaining liquid. Release ∇. Piston moves up. |
| Δ key green and flashes: | Sequence not yet finished.
Repeat points 2) – 4) to execute the remaining dilutions. |
| Δ key green constantly: | Sequence finished.
Acoustical signal: one long tone.
Sequence can be repeated immediately beginning with point 2). |

After every transfer, change tip to avoid carryover (eject tip by pressing R).

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The liquid dispensed from a Combitip is measured and displayed. Aspirating and dispensing takes place as long as the $\Delta \nabla$ keys of the dispensing grip are held down. To ensure that the end point of titration can be determined exactly, the dispensing speed is decreased by one step every time the ∇ key is pressed.

This function can only be executed with a Combitip!

Titr.Vol. parameter: Display of the volume dispensed

Programming:

Activate **command line**, select ***Titr*** function.

Activate **entry field** (see Sec. 3.2.1).

Where appropriate, enter new ***Tip*** and ***Speed*** parameters (select parameter with the "Enter" key, set values with the cursor keys). It is advisable to enter the highest speed first.

Executing with the dispensing grip:

Keys light up:

- | | |
|---------------------------------|---|
| 1) R key red : | Press R. Insert Combitip. |
| 2) Δ key green : | Filling: press Δ and hold down until enough liquid is aspirated. |
| 3) ∇ key yellow : | Liquid is dispensed as long as ∇ key is pressed. If ∇ key is pressed again, the dispensing speed is reduced by one step and remains constant when the lowest step is reached. |

If the Combitip is not sufficiently filled, aspirate liquid again. This does **not** influence the display of the volume dispensed.

Reset does **not** influence the display of the volume dispensed.

To begin a new titration (new counting of volume dispensed): press ***Mode*** and then ***Enter*** on the control unit.

5.3 *Mixing* (*Mix* function)

Mixing liquid by aspirating and dispensing with a dispensing tip.

Function can be executed with all dispensing components.

Mix.Vol. parameter: Display of the volume aspirated/dispensed to be mixed

Repeats parameter: Number of mixing strokes executed

Programming:

Activate **command line**, select **Mix** function.

Activate **entry field** (see Sec. 3.2.1).

Enter new **Tip**, where appropriate, as well as **Mix. Vol.**, **Repeats (0 – 20)**, **Speed** (select parameter with the "Enter" key, set values with cursor keys).

Executing with the dispensing grip:

Keys light up:

- 1) **R key red:** Press **R** and hold down when Pip 1000, Pip 100, Pip 10, Multi 200, Multi 50 used. Insert dispensing component and attach pipette tips, where appropriate.
- 2) **Δ key green:** Aspirating and dispensing: Hold tip in liquid, press Δ.
- 3) **Δ key green** constantly: Mixing finished.
Acoustical signal: one long tone.
Mixing can be repeated immediately beginning with point 2).

When mixing with Pip 1000, Pip 100, Pip 10, Multi 200 or Multi 50, residual liquid must be removed from the tip by executing a blow-out (as for **Pip** function). Key is yellow.

Eject the pipette tip(s) by pressing R key.

5.2.2 Pipetting with Combitips

As for **Disp** function yet liquid is aspirated for one pipetting only and is pipetted in the next step.

With Combitips, **larger** volumes can be transferred than with pipetting adapters.

Pip.Vol. parameter: Display of the volume pipetted
Repeats parameter: Repeat number of the volume pipetted

Programming, see 5.2.1.

Executing with the dispensing grip:

Keys light up:

- 1) R key **red**: Press R. Insert Combitip.
- 2) Δ key **green**: Aspirating: press Δ .
- 3) ∇ key **yellow**: Dispensing: press ∇ .
 Δ key **green** and flashes: Sequence not yet finished.
When same sample is repeatedly pipetted: repetition of points 2) and 3).
 Δ key **green** constantly: Sequence finished.
Acoustical signal: one long tone.
Sequence can be repeated immediately with point 2).

If different samples should be dispensed with the same Combitip, rinse Combitip between dispensings:

Aspirate a volume of the next sample, discard subsequently with Reset (into waste tank).

Then from Zero Position (= Reset Position, see Sec. 4.1) aspirate volume of sample and dispense.

The counting of the steps dispensed is not influenced.

Executing with the dispensing grip:

Keys light up:

- 1) R key **red**: Press R and **hold down**. Insert dispensing component. **Release R**. Attach pipette tip(s). Hold dispensing grip vertically. Immerse pipette tip 2 or 3 mm into the liquid.
- 2) Δ key **green**: Aspirating: press Δ .
- 3) ∇ key **yellow**: Dispensing:
Hold pipette tip against the inside of the vessel to be filled. Press ∇ .

Perform blow-out (remaining liquid emptied out of tip):

- 4) ∇ key **yellow**: Blow-out: hold tip against the inside of the vessel. Press ∇ and hold down. Piston moves down.
Pull tip away from the wall, release ∇ , piston moves up.
- 5) Δ key **green** and flashes: Sequence not yet finished.
Repeat points 2) – 4) to pipette the remaining number of volumes.
- 6) Δ key **green** constantly: Sequence finished.
Acoustical signal: one long tone.
Sequence can be repeated immediately beginning with point 2).

After every transfer, exchange tip to avoid carryover (eject tip by pressing R).

When pipetting viscous solutions with pipetting adapters, it is necessary to rinse the tip as when working with mechanical pipettes.

Correct rinsing:

Aspirate liquid volume into the tip. Empty tip against the inside of the vessel outside the liquid (perform blow-out). Carry out rinsing 2 or 3 times.

Pipetting small volumes ($\leq 20 \mu\text{l}$) by rinsing in a liquid, see Sec. 5.5, Example 1.

Example:

Dispensing the volume sequence
50, 100, 150, 200, 250, 300 µl, dispense each volume 3 times.

Key sequence at control unit:

After selection of **Multi Dispense** and entry of the required Combipip (e. g. 2.5 ml Combipip).

- | | Meaning: |
|--------------------------|--|
| 1) Enter: | Confirm 2.5 ml Combipip. Index number 1 is marked. |
| 2) Enter: | Activate Disp.Vol. parameter. |
| 3) Cursor keys: | Set volume 1. |
| 4) Enter: | Activate Repeats parameter. |
| 5) Cursor keys: | Set Rep. (1) = 3 |
| 6) Enter: | Activate Speed parameter. |
| 7) ▲▼cursor keys: | Set desired speed.
The set value applies to the entire sequence. |
| 8) Enter: | Index number 2 is marked, etc. |

For index number 7 for **Repeats** parameter, enter 0 (volumes 7 - 12 are skipped).

It is not possible to program a **Multi Dispense** sequence for **different** Combipips.

It is also not possible to program a different dispensing speed for every partial volume, however, while executing the sequence it is possible to re-set the speed before every dispensing (see Sec. 6.5).

Executing with the dispensing grip: as for 5.1.1

5.2 Pipetting (Pip function)

Aspirating and dispensing a liquid volume.

This function can be executed with all dispensing components.

5.2.1 Pipetting with pipetting adapters Pip 1000, Pip 100, Pip 10, Multi 200, Multi 50 (pipetting with blow-out)

- | | |
|----------------------------|--------------------------------------|
| Pip.Vol. parameter: | Display of the volume pipetted |
| Repeats parameter: | Repeat number of the volume pipetted |

Programming:

Activate **command line**, select **Pip** function.

Activate **entry field** (see Sec. 3.2.1).

Enter new **Tip**, where appropriate, as well as **Pip.Vol.**, **Repeats (0 - 99)**, **Speed** (select parameter with "Enter", set values with the cursor keys).

5.1.2 Multi Dispense = Dispensing of a sequence of different partial volumes

Disp.Vol. parameter: Display of the partial volume dispensed (shown by index)

Repeats parameter: Repeat number of the partial volume (shown by index)

With this function, you can dispense up to 12 different partial volumes in succession from one Combitip. Each partial volume can be repeated several times.

Activate **command line**, select **Disp** function.

Activate **entry field** (see Sec. 3.2.1).

With a cursor key, select **Multi Dispense**, press Enter.

An **index line** with numbers 1 – 12 appears in the entry field. Number 1 is marked (see Fig. 9), the **Disp.** parameter (= Disp.Vol.) and **Rep.** (= Repeats) are provided with index 1.


Fig. 9: Display for Multi Dispense

* = Index line. Repeats = 0 set from index number 7

! = Number of remaining dispensing steps from partial volume 1

Select new **Tip**, where appropriate, and enter a dispensing sequence of different partial volumes or modify existing sequence.

This can begin in sequence from number 1. Alternatively, the partial volume can be located on the index line with the right/left cursor keys. When an index number is marked, the parameter values **Disp.Vol.** and **Repeats** are displayed for the appropriate partial volume. A sequence already generated can thus be quickly checked and modified.

If 0 is entered for **Repeats**, this and all following steps are deleted (see Fig. 9).

5 The functions

5.1 Dispensing (*Disp* function)

Repeated dispensing of partial volumes of an amount of liquid.

This function can only be executed with Combitips!

5.1.1 Single Dispense =	Repeated dispensing of the same partial volume
Disp.Vol. parameter:	Display of the partial volume dispensed
Repeats parameter:	Repeat number of the partial volume

Programming:

Activate **command line**, select **Disp** function.

Activate **entry field** (see Sec. 3.2.1).

Set **Single Dispense**.

Select new **Tip**, where appropriate, as well as **Disp.Vol.**, **Repeats (0 – 99)**, **Speed** (select parameter with the “Enter” key, set values with the cursor keys).

Executing with the dispensing grip

Keys light up:

- | | |
|---------------------------------|---|
| 1) R key red : | Press R . Insert Combitip. |
| 2) Δ key green : | Aspirating: press Δ. |
| 3) ∇ key yellow : | Dispensing: press ∇ repeatedly. |
| Δ key green and flashes: | Combitip contains less liquid than required for 1 partial volume. Sequence not yet finished. Repeat points 2) and 3) to dispense the remaining partial volumes. |
| Δ key green constantly: | Sequence finished.
Acoustical signal: one long tone.
Sequence can be repeated immediately beginning with point 2) (see Sec. 6.1). |

Fig. 7: Inserting a dispensing component in the dispensing grip

Fig. 8: Handling the dispensing grip

4.2 Inserting the dispensing components (Fig. 7)

Press the R key of the dispensing grip. Hold down the key, if Pip 1000, Pip 100, Pip 10, Multi 200 or Multi 50 is entered as the dispensing component.

Push the barrel of the dispensing component into the slot at the lower end of the dispensing grip. Turn slightly until it clicks into position. It is not necessary to press the piston into the piston mounting.

The white adapter must be attached to the 25.0 and 50.0 ml Combitips before insertion into the dispensing grip. Adapters are included in the accessories of the device.

4.3 Handling the dispensing grip (Fig. 8)

Take the dispensing grip in your right or left hand as shown in Fig. 8. Use your thumb to press the keys.

4 The dispensing grip


Fig. 6: Dispensing grip

4.1 The control keys (Fig. 6)

R

Reset

Pressing the R (Reset) key:

The motor moves the piston mounting K to its lowest position (see Fig. 6). This is the Reset Position = Zero Position. If a dispensing component is attached, the piston is pushed completely into the dispensing component by this means.

Hold down the R key when you want to insert or remove a dispensing component: (Pip 1000, Pip 100, Pip 10, Multi 200 or Multi 50) and the name of the dispensing component is entered at the **Tip** parameter.

In this case, the piston moves back a little into the dispensing grip when the key is released (to move immediately back to the Zero Position when the R key is pressed).

Aspirating/dispensing

Depending on which of the keys have been pressed, the motor moves the piston mounting up or down.

3.2.3 Programming the control unit

After selection of a new function, the dispensing component last used and the parameter values entered are displayed in the entry field.

If a new dispensing component is set, the parameter values last entered for this component appear simultaneously.

This means that if the device is appropriately programmed, it is possible to continue working immediately with a new function/dispensing component after setting and confirmation.

The selected function can be executed as soon as the entry field is activated. The **Volume**, **Repeats** and **Speed** parameters need not be confirmed; working can continue immediately with a new value.

3.2.4 New start of a repeat count

The counting of dispensed volume begins again:

- after the volume value or the repeat number has been changed,
- after "Mode" and "Enter" have been pressed in succession.

3.3 Displaying the piston position (SUM)

The position of the piston in the dispensing component is displayed in the bottom line of the display to the right:

- as a volume value, related to the dispensing component used,
- graphically to the right highlighted on a dark background.

The volume is calculated from the home position of the appropriate function (see Sec. 6.1 and 6.2). Exception: **Manu** function with Combipip: the volume is calculated from the Zero Position (see Sec. 4.1).

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The volume to be dispensed is entered:

- with a maximum of 4 digits in μl or
- with up to 2 positions before and after the decimal point in ml.

With the $\blacktriangleleft\blacktriangleright$ cursor keys the decimal place is set. With the $\blacktriangle\blacktriangledown$ cursor keys, the figures set. If the $\blacktriangle\blacktriangledown$ cursor keys are held down, the figures scroll up to 9 and down to 0.

Each dispensing component has a minimum and maximum adjustable volume (see Sec. 2.2).

The name of the volume corresponds with the function selected, i.e. **Disp**, **Pip**, **Mix**, **Dil**, **Tit** and **Manu**.

Repeats:

- Number of repeats (repeat number) with which a volume should be dispensed, pipetted or diluted (max. 99).
- Number of mixing strokes with which liquid is aspirated/dispensed with a tip during mixing (**Mix** function, maximum number = 20).

During this function, the number of **remaining** dispensings is displayed to the right of the repeat number (see Fig. 9).

If 0 is entered as the repeat number:

- the function can be repeated as often as desired and the number of dispensings executed is displayed next to the 0 (**Single Dispense**, **Pip** and **Dil** functions, see Fig. 5),
- the appropriate step and all following steps of a dispensing sequence are not executed (**Multi Dispense** function).

Speed:

Eight aspirating/dispensing speeds are available from slow (one block) to quick (8 blocks), see Sec. 6.4.

Speed \uparrow : aspirating

Speed \downarrow : dispensing

Single Dispense/Multi Dispense (only displayed when **Disp** function selected):

With **Single Disp**, one selected partial volume can be dispensed repeatedly after filling.

With **Multi Disp**, a dispensing sequence consisting of different partial volumes can be programmed (see Sec. 5.1.2).

Strictly speaking, this parameter does not belong to the entry field, it is an extension of the command line. The **Single Dispense** or **Multi Dispense** function cannot be executed until confirmation with "Enter".

3.2 The display (Fig. 4 + 5)

The display is divided into:

a **top line**, the **command line**, in which the functions executed by EDOS can be selected,

a **middle field**, the **entry field**, in which the **parameter values** needed for dispensing are entered and displayed,

a **bottom line** in which the position of the piston in the dispensing component can be displayed on the right and a comment on the left.

3.2.1 Command line and entry field

Press "Mode":

The top line of the display, the **command line**, is activated. The line is dark, the field of the **function last executed** is highlighted (see Fig. 4).

With the ◀▶ cursor keys, a function can be selected.

Press "Enter":

The choice of function is confirmed, the **entry field** is activated.

On the bright command line, only the **field of the chosen function is dark** (see Fig. 5).

To re-activate the command line, press "Mode".

3.2.2 Entering/Displaying the parameter values in the entry field

After activation of the entry field, the Tip parameter is boxed (for the **Disp** function, press "Enter" again).

You can now execute the chosen function as defined by the set of parameter values in the entry field.

To enter new parameter values:

- select the parameter to be changed by pressing "Enter" repeatedly until it is boxed,
- set the desired parameter value with the **cursor keys**.

Tip:


The dispensing component used is set with the ▲▼ cursor keys. Press "Enter" to confirm this setting.

For the remaining parameters of the entry field, the **values last entered** now appear: the values last entered are stored automatically for each function and dispensing component. They remain in the memory even after the display has been switched off or the device disconnected from the main power supply.

By calling the appropriate dispensing component, the function can be executed immediately with the set of parameter values last used.

3 The control unit

3.1 The keypad

Press this key briefly to switch on/off the display of the control unit and the power supply to the dispensing grip.

With "Mode" and "Enter", the various fields of the display can be activated and the functions in the command line selected. With the cursor keys ▲▼◀▶, new values can be set in the entry field of the display.

When the device is switched on, the three keys of the dispensing grip light up several times (red, green, yellow).

The title of the EDOS program appears briefly in the display of the control unit, first dark on a light background and then vice versa.

The function selected before the device was switched off is then displayed. This can be executed immediately.


Fig. 4: EDOS display,
command line activated (dark)


Fig. 5: EDOS display,
entry field Disp function activated (command line highlighted)
! = Number of dispensings performed

Table 1 below lists a summary of the minimum and maximum volume which can be dispensed.

Dispensing component	Min. vol. μl	Max. vol. μl	Step size μl	Max. step no.
Pip 10	0,5	10	0.1	1
Pip 100	10	100	0.1	1
Pip 1000	100	1,000	1,0	1
Multi 50	5	50	0.1	1
Multi 200	25	200	0.2	1
0.05 ml Combitip	1	49	0.1	49
0.5 ml Combitip	5	500	0.5	100
1.25 ml Combitip	10	1,250	1.0	125
2.5 ml Combitip	20	2,500	1.0	125
5.0 ml Combitip	40	5,000	1.0	125
12.5 ml Combitip	100	12,500	10.0	125
25.0 ml Combitip	200	25,000	10.0	125
50.0 ml Combitip	400	50,000	10.0	125

Min. Vol. = Minimum volume per dispensing step.

Max. Vol. = Maximum volume per dispensing step.

Step size μl = Smallest volume difference adjustable.

With Multi 200 and 0.5 ml Combitip, 0.1 μl steps can be set and dispensed.

After the first keypress, the values rounded up to the corresponding step sizes are displayed.


Fig. 3: Dispensing components

- a) 0.05, 0.5, 1.25, 2.5, 5.0, 12.5, 25.0 and 50.0 ml Combitips
- b) Pip 1000, Pip 100, Pip 10
- c) Multi 200 und Multi 50
- d) 1,000, 100 und 10 µl pipette tips

2.2 The dispensing components (Fig. 3)

The liquid is filled into dispensing components which are inserted into the dispensing grip:

Eppendorf Combitips,

all sizes (0.05 ml to 50.0 ml).

The size of the Combitip indicates the maximum volume which can be aspirated. Liquid is aspirated **directly into the Combitips**.

A 10 µl pipette tip is attached to the 0.05 ml Combitip. Liquid is **only** aspirated **into the pipette tip**.

Pip 1000, Pip 100, Pip 10

These are pipetting adapters (lower pipette housings) to which 1,000 µl, 100 µl and 10 µl pipette tips are attached. Liquid is **only** aspirated **into the pipette tips**.

Multi 200 and Multi 50 (in preparation)

Pipetting adapters for multi-channel dispensing (lower housings of 8-channel pipettes) for 100 µl pipette tips. With these adapters, a maximum of 200 or 50 µl can be pipetted per channel. Liquid is **only** aspirated **into the pipette tips**.

In many cases, Combitips enable dispensing in free jet with exact drip separation. To attain optimal values, it may be necessary to attach a 100 µl pipette tip to the Combitip and hold the tip against the inside of the vessel (see Sec. 11).

Components of the device


Fig. 1: EDOS

2.1 EDOS (Fig. 1 + 2)

The electronic dispensing system EDOS consists of:

- the dispensing grip (1) with keys for initiating aspirating and dispensing,
- the control unit (2) with keys for selecting the functions and entering the volumes and parameters,
- power unit for the power supply (3).

Before plugging in the device, compare your local power supply with the electrical requirements listed on the identification plate on the power unit.

Connect power unit, control unit and dispensing grip as shown in Fig. 1. Begin by connecting the dispensing grip and the control unit. To avoid unintentional loosening of the connections, secure the plugs with the screws and the screw sleeve.

The rest for the dispensing grip (4) is placed onto the right side of the control unit (Fig. 2).

The brightness of the display can be set with the knurled screw on the right side of the housing (5).


Fig. 2: EDOS control unit

1 Introduction

Dear Customer,

The electronic dispensing system EDOS from Eppendorf is an all-purpose device for dispensing, pipetting, mixing, measuring and titrating liquids:

- **Dispensing** liquid from a filled dispensing tip in whatever partial volumes desired: **Disp** function
- **Pipetting** with ejectable tips (as with manual pipettes): **Pip** function
- **Mixing** liquids by successive aspirating and dispensing with a dispensing tip: **Mix** function
- **Titration** by exactly measuring a liquid volume dispensed: **Tit** function
- **Dilution series** executable in a programmable order: **Dil** function
- **Aspirating** and **dispensing** of liquid without preselected volume: **Manu** function.

The functions of the device are initiated by hand as for manual pipettes but without applying force. The intelligent software of the device supports their execution.

Using the well-known plastic dispensing tips from Eppendorf (Combitips, pipette tips), a variety of samples and liquids can be dispensed in succession simply by exchanging the tip.


When used correctly with Eppendorf dispensing tips, EDOS permits liquid measurement with high accuracy and precision.

For applications which demand especially high purity, we recommend biopur pipette tips and Combitips which are colorless, sterile, pyrogen-free, DNA-free, RNase-free and ATP-free.


Due to these features, EDOS 5221 is suitable for use in biochemical, molecular biological, clinical and industrial laboratories.

Contents

1	Introduction	2
2	Components of the device	3
2.1	EDOS	3
2.2	The dispensing components	4-5
3	The control unit	6
3.1	The keypad	6
3.2	The display	7
3.2.1	Command line and entry field	7
3.2.2	Entering/Displaying the parameter values in the entry field	7-8
3.2.3	Programming the control unit	9
3.2.4	New start of a repeat count	9
3.3	Displaying the piston position (SUM)	9
4	The dispensing grip	10
4.1	The control keys	10
4.2	Inserting the dispensing components	11
4.3	Handling the dispensing grip	11
5	The functions	12
5.1	Dispensing (Disp function)	12
5.1.1	Single Dispense	12
5.1.2	Multi Dispense	13-14
5.2	Pipetting (Pip function)	14
5.2.1	Pipetting with pipetting adapters	14-15
5.2.2	Pipetting with Combipips	16
5.3	Mixing (Mix function)	17
5.4	Titration (Titr function)	18
5.5	Diluting (Dil function)	19-20
5.6	Aspirating and dispensing any amounts of liquid (Manu function, manual mode)	21
6	Useful hints	22
6.1	Liquid transfer with a Combipip	22-23
6.2	Liquid transfer with pipetting adapters	23
6.3	Controlling the piston movement, Reset	24
6.3.1	Stopping piston movement	24
6.3.2	Piston movement with the Titr and Manu functions	24
6.3.3	Reset function	24
6.4	Speed of piston movement ("Speed" parameter)	25
6.5	Entering parameter values during execution of a function, new start of a repeat counting	25
6.6	User prompting	26
6.6.1	Optical signals	26
6.6.2	Acoustical signals	27
6.6.3	Comments	27
7	Storing and selecting programs (Pgm function)	28
7.1	Storing a set of parameter values	28
7.2	Selecting a set of parameter values	29
8	Options (Opt function)	30
9	Troubleshooting	31-32
10	Maintenance, cleaning, spare parts	33-34
11	Technical data	35-36
12	Ordering information	37-38

EDOS® 5221
Instruction Manual

Eppendorf –
Netheler – Hinz GmbH
D-22331 Hamburg
Tel.: (0 40) 5 38 01-0
Teletex: 40 30 61 = EGHAM
Fax: (0 40) 5 38 01 556

eppendorf

eppendorf® is a registered trademark